

The A-Z of *Azlon*[®]

REUSABLE PLASTIC LABWARE




Azlon®

For more than half a century, the Azlon® name has been synonymous with high quality affordable plastic laboratory ware.

Available in a variety of polymers, Azlon® reusable plastics cover a broad range of applications in the modern laboratory, offering more than 700 products, including bottles, wash bottles, measuring cylinders, jugs and beakers.

Azlon® products are manufactured to the highest quality standards utilising a range of plastic moulding techniques including blow and injection moulding.

Leadership in polymer science and many important innovations have combined to make Azlon® one of the most widely known and respected brands in durable plastic labware.



Azlon® key products

Container bottles

Azlon® leak-proof container bottles provide a safe and secure storage option that you can trust!

The linerless design prevents contamination whilst a smooth inner neck and base contours allow easy pouring and access for cleaning.

Beakers

Azlon® beakers are available in two styles: classic, straight-sided for ultimate stability and tapered, stackable for efficient use of storage space.

Volumetrics

Perform precise measuring with the Azlon® volumetrics range! From volumetric flasks to measuring cylinders all made from virtually unbreakable polymers, you can enjoy the same accuracy as glassware without the safety hazard.

Wash bottles

The market leading Azlon® multilingual wash bottle range is manufactured from high quality LDPE, allowing maximum flexibility whilst maintaining excellent chemical resistance.

Patented Driplock® vapour venting technology aids a precise and continuous stream of liquid, even in the case of volatile solvents.

GHS compliant pre-labelling with chemical identification and handling information ensures the safe and accurate dispensation of liquids.



Azlon® wash bottles

It's no surprise that Azlon® wash bottles are the market leaders!

As well as having the biggest selection of pre-labelled bottles, the Azlon® custom print service makes designs in every language possible and patented Driplock® vapour-venting technology ensures safe and accurate dispensing of chemicals.

- **Multilingual:** English, French, German, Scandinavian, Polish, Arabic, Spanish, Russian.
- **Driplock® Vapour Venting Technology:** Eliminates loss of volatile solvents by preventing evaporation and stops pressure build ups, putting an end to hazardous, costly drips and ensuring an accurate and continuous dispensing stream.
- **Custom Print Service:** Design your own artwork for printing onto the Azlon® wash bottle of your choice.



Trackable with complete labelling facility ensuring full traceability

DIY write-on-panel for reusability

Slim design for ease of handling

Key Facts

Competition	<ul style="list-style-type: none"> • Bel Art, Kartel, Vitlab
Advantages over Competition	<ul style="list-style-type: none"> • Available in more languages • Custom printing • Cost-saving, patented DripLock system • Broadest range on the market
Price position	<ul style="list-style-type: none"> • Low
Application/Industry use	<ul style="list-style-type: none"> • Industrial Cleaning • Scientific & Chemical • Clinical & hospitals

Attention! Achtung! انتباه



Azlon® volumetrics

The Azlon® volumetrics range has all the essentials needed for measuring in the lab, from conical and volumetric flasks to measuring cylinders and burettes – and all with the same accuracy as glass but without the hazard!


- Class A & Class B tolerances for accuracy
- PMP material for glass-like clarity
- Easy-to-read markings for enhanced usability



Key Facts

Competition	<ul style="list-style-type: none"> • Kartel, Vitlab, Brand
Advantages over Competition	<ul style="list-style-type: none"> • Draining hole in measuring cylinder • No-slip, anti-suction base • Moulded and over-printed graduations
Price position	<ul style="list-style-type: none"> • Mid-range
Application/Industry use	<ul style="list-style-type: none"> • Scientific & Chemical • Clinical & hospitals • Petrochemical

made to measure



Azlon® beakers

Style	Size	Polypropylene	PMP	Printed	Moulded	Autoclavable
Classic – straight-sided	10ml to 2litre	✓	✓	✓	✗	✓
Tapered – stackable	25ml to 5litre	✓	✗	✓	✓	✓
Jugs	50ml to 5litre	✓	✗	✓	✓	✓

Perform accurate pouring with Azlon® beakers!
 Available in straight-sided, stackable and jug versions,
 the Azlon® beaker range gives you the perfect long term
 solution for stirring, mixing and measuring.



Key Facts

Competition	<ul style="list-style-type: none"> • Kartell, Vitlab, Brand
Advantages over Competition	<ul style="list-style-type: none"> • Unique mouth profile for maximum grip • Available in multipacks • Maximum access to base for wide equipment
Price position	<ul style="list-style-type: none"> • Mid
Application/Industry use	<ul style="list-style-type: none"> • Food & Beverage • Scientific & Chemical • Clinical & hospitals

**the must-have
tool for any lab**

precise non-drip
pouring spout

Unique mouth
profile for
maximum grip

Crisp, blue,
easy-to-read
graduations



Azlon® storage bottles

Trust Azlon® with your storage needs! Azlon® storage bottles are available in many shapes and sizes and with a leakproof, linerless design that ensures a secure closure. Manufactured using the highest quality resins, the various polymer options available ensure you will find the perfect solution for your application.

- **Leakproof*:** High quality moulded threads and seals ensure secure closure.
- **Linerless:** Eliminates contamination from loose cap insert liners.
- **Base mouldings:** Include recycling code, bottle volume and registration ramp.



*Leakproof when specified in product description.

Key Facts

Competition	<ul style="list-style-type: none">• Nalgene, Burkle
Advantages over Competition	<ul style="list-style-type: none">• High quality whilst retaining affordability• Broad range of shapes sizes and polymers
Price position	<ul style="list-style-type: none">• Mid
Application/Industry use	<ul style="list-style-type: none">• Food & Beverage• Scientific & Chemical• Clinical & hospitals

**secure & safe ...
storage you
can trust**

Wide neck for
easy filling

Pronounced collar
for anchoring
shrinkwrap

Narrow-neck for
accurate pouring



Write-on-panel for
easy labelling



The A-Z of
Azlon
REUSABLE PLASTIC LABWARE

Get in touch if you have any questions
and please remember...

We want to hear from you!!!

General email: azlon@scilabware.com
Product Manager email: haneen.al-azzawi@scilabware.com

SciLabware Limited

Unit 4, Riverside 2
Campbell Road
Stoke-on-Trent
Staffordshire
ST4 4RJ
United Kingdom
Tel: +44 (0)1782 444406
Fax: +44 (0)1782 940436

Registered in England and Wales No. 7723574
BS EN ISO 9001:2008 registered company



SciLabware